



### Efficiently Mobile

The SideCar is GlobalMed's enhancement to the Polycom Practitioner Cart for mobile video conferencing. It stores and connects GlobalMed medical peripherals in keeping with the design of the Practitioner Cart base.

- The SideCar, when added to the Practitioner Cart, integrates medical peripherals and respective cables to reduce clutter and possible cord damage.
- Preserves the sleek look of the Practitioner Cart in an easy-to-clean configuration.
- SideCars work as splash guards to protect interior peripherals from contamination.

### Fully Functioning Telemedicine Platform

- One SideCar can be attached to each side of the Practitioner Cart's base.
- SideCars house GlobalMed's TotalENT™ otoscope and a PC containing GMGs CapSure® Image Automation Software.
- GlobalMed's TotalExam™ camera and retractable cable reel can be mounted on the exterior of the SideCar.
- Both SideCars offer doors that provide accessibility to the devices housed inside.



### More Uses While Maintaining a Small Footprint

The SideCar adds greater functionality without interfering with height adjustment. Users benefit from the added medical peripherals without impeding workflow.

- Each SideCar is constructed of steel for structural durability.
- Incorporating two SideCars to a dual-monitor cart does not increase overall footprint.
- Rear of SideCars are enclosed to prevent accidental disconnection of peripherals.

Works with your PC or Mac and on most popular video codecs including:





**SideCar Dimensions**

Width: 4.5" per SideCar  
Depth: 18.25"  
Height: 17"

**TotalExam™ 2 Camera S-Video**

Lightweight.  
Up-Close Imagery.  
Ergonomically Designed.  
LED Light Carousel.

**TotalENT™ Combination Camera and Light Source**

Crystal Clear German Optics.  
Crisp Image.  
Easy Telemedicine Integration.

**High-End PC**

Small form factor  
HD graphics and output  
Incorporates with clinical IT systems.

**CapSure®**

Patient Safety.  
PACS Integration.  
DICOM / HL7 Workflow.



Works with your PC or Mac and on most popular video codecs including:

