



WallDoc® is a telemedicine solution designed for areas with space limitations. WallDoc mounts flat on the wall and occupies the same amount of space as a flat panel monitor, yet provides the functionality of a fully-integrated telemedicine solution. From a small clinic, corporate health center, or school nurse's office to an oil rig, WallDoc is the comprehensive, stationary telemedicine platform for those with space challenges.

Efficient and Intelligent

- Simple to launch and use
- Customizable to your workflow
- Compatible with connected instruments
- Integrated software and hardware

Secure and Compliant

- Integrated security features
- Prevents unauthorized usage

Streamlined Clinical Aesthetic

- Sleek, ergonomic design
- Adjusts to standing or seated positions
- Stores connected medical devices and exam cameras
- eNcounter® Software customizable to your brand

For more information on WallDoc and accompanying peripheral devices, call **+1 (877) 959-4713** or visit **GlobalMed.com**.



Component Specification	Description
Dimensions	26" H, 26" W, 16" D minimum 32" H, 38" W, 22" D maximum with main camera, slide out device compartment and fully extended keyboard tray
Computer	Integrated Dell computer
Operating System	Windows® 7 or higher
Display	27" 1080p integrated monitor
Surface	Anti-microbial casing
Keyboard	Wireless full-size keyboard and mouse
Compartments	Secure integrated device storage, Security power lock and Slide out device compartment tray
Ports	10 USB 3.0 ports
Audio	Integrated Speaker
Video Conferencing	Supports hardware and software video conferencing
Device Integration	GlobalMed approved medical devices (call for validated list)
Movement	Vertical: 29.5" - 42" Standing height Extension: 22" max from wall