



TOTALEXAM® Backpack

The GlobalMed® Transportable Exam® Backpack (TEB) brings a new dimension to providing telemedicine to the patient in remote environments.

This lightweight, USB-powered design provides a simple and organized solution, configured for in-home care or remote provider consultations with our suite of clinical devices including:

- TotalVitals® Health Monitor
- TotalExam® Lite, Camera & Otoscope
- USB Stethoscope
- TotalECG®
- Ultrasound

Our devices allow a wide array of clinical data collection, from vitals signs including blood pressure, heart rate and temperature, to heart and lung sounds. Even in the most rural areas or isolated locations, a clinician can gather video evidence for an ENT consultation, and even capture ultrasound images, expanding care to patients anytime, anywhere.



+1 (480) 922-0044
globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260
United States



TOTALEXAM® Backpack



Items in the TotalExam® Backpack may be ordered separately. The complete kit includes:

TotalExam® Lite Camera | Qty 1
TotalExam® Lite Oto | Qty 1
TotalExam Polarizing Derm Hood | Qty 1
Tongue Depressor Adapter | Qty 1
USB 2.0 Cable, 2.5 Meters | Qty 2
TotalVitals® | Qty 1
ClearSteth® Stethoscope | Qty 1
TotalECG® | Qty 1
Blood Pressure Cuffs | Qty 3
Ultrasound | Qty 1
Headphones | Qty 1
One Touch Collars | Qty 5
6" Adult Tongue Depressor | Qty 5
Adult Specula, 4.25mm | Qty 5
Child Specula, 3.0mm | Qty 5

FEATURES

- Lightweight and compact bag design
- Integrated 4 port USB hub
- Dedicated laptop pocket up to 15.6"
- Custom foam insert to support peripheral devices.
- Additional storage pockets for peripheral accessories