



CLINICALACCESS® Station Lite

ClinicalAccess® Station Lite (CAS Lite) frees you from the space concerns of larger form factors. Customize your station to accommodate your examination needs using its innovative, modular design. With its slim footprint, built-in WiFi and Ethernet, CAS Lite is the perfect mobile telemedicine solution to seamlessly fit into small clinical spaces.

Efficient and Intelligent

- Simple to launch and use
- Customizable to your workflow
- Integrated software and hardware

Secure and Compliant

- Integrated security features
- Locking compartment option, prevents unauthorized usage

Sleek and Mobile

- Sleek, small form factor design
- Small footprint and five-wheel base with 5" casters for stability, ease of egress, and turning
- Highly maneuverable in confined spaces
- Spring-assisted adjustable ergonomics; worksurface movement is 33-44"


+1 (480) 922-0044
globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260
United States





CLINICALACCESS® Station Lite



Dimensions	20"W x 24"L x 53"H (lowest position, excluding camera)
AC Power	90V to 220V, 50Hz to 60Hz Voltage: Nominal 120V Input Voltage Range: 90VAC- 148VAC Frequency: 60Hz ± 3 Hz Plug Type: NEMA 5-15P
Battery Power	Dual battery Range of 6-8 Hours (varies by usage) Sealed Lead-Acid (SLA) VA 1000VA, Watts 600W, Size 12V/9AH Typical Recharge Time: 8 hours
Outlets	Ethernet and AC power, located on back of integrated dust cover
Computer	Integrated Dell® or like model
Ports	USB hub (4 ports) 5-port switch option within base
Video Conferencing	Compatible with most major hardware and software codecs including: eNcounter®View, Cisco®, Vidyo®, Avaya®, and others. Contact the GlobalMed sales team to learn more.
Display	27" LG Monitor 22" LG Monitor Dual 22"/ 16" Monitor with swing arm
Camera	PTZ camera, standard 12x optical zoom
Sound	Dell® Sound Bar 17.4"W x 2.4"D x 1.4"H Speaker Type: Active 1 x sound bar 5 Watt - 90 - 20000 Hz
Table Top	19 1/2"W x 20" L; Table top height range is 33" to 44"
Security	Bin storage in tabletop, with combination lock
Cables	10' Ethernet cable reel, 8' power cable reel