



NEXTEXAM™ Mobile Telemedicine Station

Clean, Compact & Organized

GlobalMed, the first company to integrate video conferencing equipment and medical peripherals into one portable system takes the platform to the next level by housing all equipment in one highly organized and secure solution. The NextExam™ modernizes the look by enclosing all video conferencing equipment and respective cables out of site but still within reach of the clinician.

The NextExam™ can be equipped with your choice of medical devices to be contained within the peripheral bins. The system is available with user-selected hardware or software codec for compatibility with legacy systems.




+1 (480) 922-0044
globalmed.com

15023 N 73rd Street
Scottsdale, AZ 85260
United States



NEXTEXAM™ Mobile Telemedicine Station



Dimensions	74"x31" Base x 40" Tabletop
Mobility	6 wheels with locking castors
Security	2 keyed locks
Power	Uses 110V/120VAC 50/60 Hz @5.1A 2 self-powered USB 3.0 Ports Battery will run fully loaded for a minimum of 1 hour
Communications	Speaker and Echo-free microphone Wifi Ethernet
Operating System	Windows® 10 Professional
PC Options	HP ZBook Laptop
Display	27-inch diagonal LED-monitor FULL HD (1920x1080)
Processor	Intel® Core™ i7 9750H
Memory	4 DIMM slots: up to 32GB 1600 MHz
Internal Storage	750GB 7200 rpm SMART SATA II HDD
Wireless Support	Intel® Centrino® 802.11a/b/g/n, HP Integrated Module with Bluetooth® v2.1+ EDR Wireless Technology
Input Device	Full-sized spill-resistant keyboard, touchpad with scroll zone and gestures support
Camera	PTZ camera, standard 12x optical zoom
Designated Power Supply	TrippLite HC350SR Input 110VAC nominal 60Hz Output 110VAC nominal Tested, certified and compliant to UL 60601-1, IEC 60601-1-2 3rd edition and IEC 60601-1-1 3rd Edition

MODEL NUMBER GMR12150000