



# THE WELCH ALLYN® DIAGNOSTIC CARDIO SUITE™

Embedded into GlobalMed's eNcounter® virtual care platform

The Welch Allyn Diagnostic Cardiology Suite now includes ECG and spirometry capabilities, helping you shorten the path to diagnosis for improved patient outcomes. Simplify clinical workflows, reduce manual data entry and minimize complexity with one device that can perform two reimbursable tests. Our flexible platform allows you to adapt the Diagnostic Cardiology Suite to meet your clinical needs—whether you want to perform both exam types, or prefer a standalone ECG or spirometry solution. We're combining operational experience with clinical excellence to bring you one solution for two critical exams—because simple is smart.



## Features

- Wireless acquisition with the Welch Allyn Wireless Acquisition Module
- ECG applications launch directly from the EMR
- Captures best 10 seconds of ECG data from memory, which helps reduce need for repeat ECGs
- Utilize one device for spirometry and ECG exams for a common workflow and setup
- Supports DICOM® and HL7® communication
- Securely control user logins, assign user roles and access audit trails
- ECG complies with ACC, AHA and HRS recommendations for adult/pediatric patients
- Instantly verify spirometry exam quality and test performance
- Streamlines workflow by capturing and saving an ECG in as few as two clicks
- Easily access prior spirometry exams to monitor treatment plans over time
- VERITAS® ECG algorithm limits filtering to help show critical waveform data
- Repeatability criteria provides feedback and opportunities for exam coaching
- Save PDFs to any network folder
- Help perform spirometry exams correctly with an incentive screen

+1 (480) 922-0044  
globalmed.com

15023 N 73<sup>rd</sup> Street  
Scottsdale, AZ 85260  
United States

MODEL NUMBER GMD50500004